


Anti-Nephrin Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O60500

Item Number: CM0004472



Introduction

This rabbit monoclonal antibody targets human Nephrin (SWISS: O60500), validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying podocyte slit diaphragm and kidney filtration.

[Learn More](#)

Parameter	Value
Product Name	Nephrin Rabbit mAb
Remarks/Alias	CNF; NPHN; nephrin; Nephrin
Species	Human
Gene ID (Human)	4868
Gene ID	4868
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1019-1241 of human Nephrin (NP_004637.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60500
Protein Weight	135kDa
Shipping	Ice bag