

Anti-Nesprin 1 Monoclonal Antibody For Wb, If/Icc, And Elisa - Q8NF91

Item Number: CM0003783



Introduction

Rabbit monoclonal antibody against human Nesprin 1 (SYNE1), validated for WB, IF/ICC, and ELISA, cross-reacting with human and mouse. Ideal for studying LINC complex and nuclear envelope dynamics.

[Learn More](#)

Parameter	Value
Product Name	Nesprin 1 Rabbit mAb
Remarks/Alias	8B; AMC3; AMCM; CPG2; ARCA1; EDMD4; KASH1; MYNE1; Nesp1; SCAR8; C6orf98; dj45H2.2; Nesprin 1
Species	Human
GeneID (Human)	23345
GeneID	23345
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 8200-8400 of human Nesprin 1 (Q8NF91).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NF91
Protein Weight	1011kDa/112kDa
Shipping	Ice bag