

Anti-Nesprin 1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Nf91

Item Number: CM0021105



Introduction

High-quality rabbit polyclonal antibody targeting human Nesprin 1 (SYNE1), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for cytoskeleton and nuclear envelope research.

[Learn More](#)

Parameter	Value
Product Name	Nesprin 1 Rabbit pAb
Remarks/Alias	8B; AMC3; AMCM; CPG2; ARCA1; EDMD4; KASH1; MYNE1; Nesp1; SCAR8; C6orf98; dj45H2.2; Nesprin 1
Species	Human
GeneID	23345
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 8220-8410 of human Nesprin 1 (NP_892006.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NF91
Protein Weight	1011kDa/112kDa
Shipping	Ice bag