

Anti-Nesprin3 (Kash3) Monoclonal Antibody For Wb, Elisa - Q6Zmz3

Item Number: CM0002646



Introduction

Rabbit monoclonal antibody targeting human Nesprin-3 (KASH3), validated for WB and ELISA with reactivity in human and mouse. Essential for LINC complex-mediated nuclear-cytoskeletal connections and directional cell migration.

[Learn More](#)

Parameter	Value
Product Name	Nesprin3 Rabbit mAb
Remarks/Alias	KASH3; NET53; Nesp3; C14orf49; C14orf139; LINC00341; NCRNA00341; Nesprin3
Species	Human
Gene ID (Human)	161176
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 750-935 of human Nesprin3 (NP_689805.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6ZMZ3
Protein Weight	112kDa
Shipping	Ice bag