

Anti-Nestin Polyclonal Antibody For Wb - P48681

Item Number: CM0016893

Introduction

Rabbit polyclonal antibody against human Nestin (NES) for WB, IHC-P, IF/ICC, ELISA. Suitable for neural stem cell and developmental research. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	Nestin Rabbit pAb
Remarks/Alias	Nbla00170; Nestin
Species	Human
GeneID (Human)	10763
GeneID	10763
Immunogen	Recombinant protein / a synthetic peptide corresponding to a sequence within amino acids 1392-1621 of human NES (NP_006608.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48681
Protein Weight	177 kDa
Shipping	Ice bag