

Anti-Nestin Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P48681

Item Number: CM0005104

Introduction

Rabbit monoclonal antibody against human Nestin (177kDa), validated in WB, IF/ICC, IF-P, IHC-P, and ELISA. Suitable for neural stem cell and neurodevelopment research.

[Learn More](#)

Parameter	Value
Product Name	Nestin Rabbit mAb
Remarks/Alias	Nbla00170
Species	Human
Gene ID (Human)	10763
Gene ID	10763
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1392-1621 of human Nestin (NP_006608.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P48681
Protein Weight	177kDa
Shipping	Ice bag