


Anti-Net1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q7Z628

Item Number: CM0014676



Introduction

Rabbit polyclonal antibody targeting human NET1 (RhoGEF proto-oncogene). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse. Immunogen: recombinant fusion protein (aa 1-300).

[Learn More](#)

Parameter	Value
Product Name	NET1 Rabbit pAb
Remarks/Alias	NET1A; ARHGEF8; NET1
Species	Human
GeneID (human)	10276
GeneID	10276
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human NET1 (NP_001040625.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z628
Protein Weight	68kDa
Shipping	Ice bag