

Anti-Neu1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99519

Item Number: CM0028831

Introduction

Rabbit polyclonal antibody against human NEU1 (Sialidase-1), validated for Western blot, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for lysosomal storage disease and glycobiology research.

[Learn More](#)

Parameter	Value
Product Name	NEU1 Rabbit pAb
Remarks/Alias	NEU; NANH; SIAL1; NEU1
Species	Human
GeneID (Human)	4758
GeneID	4758
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-415 of human NEU1 (NP_000425.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99519
Protein Weight	45kDa
Shipping	Ice bag