


Anti-Neu2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y3R4

Item Number: CM0030446



Introduction

High-quality rabbit polyclonal antibody targeting human NEU2 (Sialidase-2) for Western blot and ELISA. Suitable for glycolipid metabolism research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	NEU2 Rabbit pAb
Remarks/Alias	SIAL2; NEU2
Species	Human
Gene ID (Human)	4759
Gene ID	4759
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 171-380 of human NEU2 (NP_005374.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y3R4
Protein Weight	42kDa
Shipping	Ice bag