

Anti-Neu3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uq49

Item Number: CM0017220

Introduction

NEU3 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Recognizes human and mouse sialidase-3, a 48 kDa membrane sialidase involved in ganglioside catabolism and EGFR signaling.

[Learn More](#)

Parameter	Value
Product Name	NEU3 Rabbit pAb
Remarks/Alias	SIAL3; NEU3
Species	Human
Gene ID (Human)	10825
Gene ID	10825
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 222-461 of human NEU3 (NP_006647.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UQ49
Protein Weight	48kDa
Shipping	Ice bag