

Anti-Neun (Fox-3) Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - A6Nfn3

Item Number: CM0014026

Introduction

Rabbit polyclonal antibody targeting human NeuN (RBFOX3), validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for neuronal marker studies and splice variant research.

[Learn More](#)

Parameter	Value
Product Name	NeuN Rabbit pAb
Remarks/Alias	FOX3; NEUN; FOX-3; HRNBP3; NeuN
Species	Human
Gene ID (Human)	146713
Gene ID	146713
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of NeuN (NP_001076044.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	A6NFN3
Protein Weight	34kDa
Shipping	Ice bag