

Anti-Neun Monoclonal Antibody For Wb, Ihc, If - A6Nfn3

Item Number: CM0003757

Introduction

Recombinant rabbit monoclonal anti-NeuN antibody suitable for Western blot, immunohistochemistry, immunofluorescence, and ELISA detection of neuronal nuclei in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	NeuN Rabbit mAb
Remarks/Alias	FOX3; NEUN; FOX-3; HRNBP3; NeuN
Species	Human
GeneID	146713
Immunogen	A synthetic peptide corresponding to amino acids 1-100 of human NeuN (A6NFN3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IF-P, IHC-P, mIHC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	A6NFN3
Protein Weight	34kDa
Shipping	Ice bag