

Anti-Neuraminidase/Neu1 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q99519

Item Number: CM0008043

Introduction

Rabbit monoclonal antibody against human Neuraminidase/NEU1 (Sialidase-1), validated for WB, IF-P, ELISA. Cross-reacts with human and mouse. Ideal for lysosomal sialidase research and glycoprotein metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	Neuraminidase/NEU1 Rabbit mAb
Remarks/Alias	NEU; NANH; SIAL1
Species	Human
GeneID (Human)	4758
GeneID	4758
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 48-415 of human Neuraminidase/NEU1 (NP_000425.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99519
Protein Weight	45kDa
Shipping	Ice bag