

Anti-Neuregulin-3 (Nrg3) Monoclonal Antibody For Wb, Elisa - P56975

Item Number: CM0003018

Introduction

Rabbit monoclonal antibody against human Neuregulin-3 (NRG3), validated for WB and ELISA. Suitable for detecting NRG3 in human, mouse, and rat samples. A tool for investigating NRG3-ERBB4 signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Neuregulin-3 (NRG3) Rabbit mAb
Remarks/Alias	HRG3; pro-NRG3; Neuregulin-3 (NRG3)
Species	Human
GeneID (Human)	10718
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-720 of human Neuregulin-3 (NRG3) (P56975)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56975
Protein Weight	78kDa
Shipping	Ice bag