

Anti-Neurod1 Rabbit Polyclonal Antibody For Ihc-P, Elisa - Q13562

Item Number: CM0023624

Introduction

NEUROD1 Rabbit Polyclonal Antibody for IHC-P and ELISA. Recognizes human and mouse NEUROD1, a bHLH transcription factor critical for neurogenesis and pancreatic islet cell differentiation. Useful in developmental biology and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	NEUROD1 Rabbit pAb
Remarks/Alias	T2D; BETA2; BHF-1; MODY6; NEUROD; bHLHa3; NEUROD1
Species	Human
Gene ID (Human)	4760
Gene ID	4760
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-356 of human NEUROD1 (NP_002491.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	Q13562
Protein Weight	40kDa
Shipping	Ice bag