

Anti-Neurod6 Polyclonal Antibody For Wb, Elisa - Q96Nk8

Item Number: CM0021982

Introduction

Rabbit polyclonal antibody against human NEUROD6 (Nex1, Atoh2, MATH2). Validated for WB and ELISA with cross-reactivity to mouse and rat. Suitable for neurodevelopmental research and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	NEUROD6 Rabbit pAb
Remarks/Alias	Nex1; Atoh2; MATH2; NEX1M; Math-2; bHLHa2; NEUROD6
Species	Human
GeneID (Human)	63974
GeneID	63974
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 188-337 of human NEUROD6 (NP_073565.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96NK8
Protein Weight	39kDa
Shipping	Ice bag