

Anti-Neurofilament H Monoclonal Antibody For Wb, Ip, Elisa - P12036

Item Number: CM0002058

Introduction

Rabbit monoclonal antibody targeting human Neurofilament H (NFH). Validated for WB, IP, ELISA; cross-reacts with mouse and rat. A key intermediate filament protein for neuronal caliber and axon stability; ideal for ALS and CMT research.

[Learn More](#)

Parameter	Value
Product Name	Neurofilament H Rabbit mAb
Remarks/Alias	NFH; CMT2CC; Neurofilament H
Species	Human
GeneID (Human)	4744
GeneID	4744
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12036
Protein Weight	112kDa
Shipping	Ice bag