

Anti-Neurofilament H Polyclonal Antibody For Wb, Ip, If, Elisa - P12036

Item Number: CM0026391

Introduction

Rabbit polyclonal antibody against Neurofilament H (NEFH), validated for WB, IP, IF, and ELISA. Reacts with human, mouse, and rat. Suitable for neuronal cytoskeleton and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	Neurofilament H Rabbit pAb
Remarks/Alias	NFH; CMT2CC; Neurofilament H
Species	Human
Gene ID (Human)	4744
Gene ID	4744
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12036
Protein Weight	112kDa
Shipping	Ice bag