

Anti-Neurofilament L Monoclonal Antibody For Wb, If, Ihc, Elisa - P07196

Item Number: CM0003522

Introduction

High-quality Neurofilament L Rabbit monoclonal antibody validated for WB, IF, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for neuronal cytoskeleton research and neurodegenerative disease studies.

[Learn More](#)

Parameter	Value
Product Name	Neurofilament L Rabbit mAb
Remarks/Alias	NFL; NF-L; NF68; CMT1F; CMT2E; CMTDIG; PPP1R110; Neurofilament L
Species	Human
GeneID (Human)	4747
GeneID	4747
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-250 of human Neurofilament L (NP_006149.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07196
Protein Weight	62 kDa
Shipping	Ice bag