

Anti-Neurofilament L Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P07196

Item Number: CM0015436

Introduction

Rabbit polyclonal antibody targeting human neurofilament light polypeptide (NEFL). Validated for WB, IF-P, IHC-P, ELISA with cross-reactivity to mouse and rat. For neuronal cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	Neurofilament L Rabbit pAb
Remarks/Alias	NFL; NF-L; NF68; CMT1F; CMT2E; CMTDIG; PPP1R110; Neurofilament L
Species	Human
Gene ID (Human)	4747
Gene ID	4747
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-250 of human Neurofilament L (NP_006149.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07196
Protein Weight	62 kDa
Shipping	Ice bag