

Anti-Neurofilament M Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P07197

Item Number: CM0022780

Introduction

Rabbit polyclonal antibody against Neurofilament M (NFM/NEF3), validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, and rat. Suitable for neuronal cytoskeleton and intermediate filament research.

[Learn More](#)

Parameter	Value
Product Name	Neurofilament M Rabbit pAb
Remarks/Alias	NFM; NEF3; NF-M; Neurofilament M
Species	Human
Gene ID (Human)	4741
Gene ID	4741
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-260 of human Neurofilament M (NP_005373.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07197
Protein Weight	102kDa
Shipping	Ice bag