

Anti-Neurog3 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y4Z2

Item Number: CM0030915

Introduction

Rabbit polyclonal antibody against human NEUROG3, suitable for WB, IHC-P, and ELISA. Reacts with human, mouse, rat. Ideal for studying enteroendocrine differentiation and pancreatic endocrine cell specification.

[Learn More](#)

Parameter	Value
Product Name	NEUROG3 Rabbit pAb
Remarks/Alias	ngn3; Atoh5; NGN-3; Math4B; bHLHa7; NEUROG3
Species	Human
Gene ID	50674
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human NEUROG3 (NP_066279.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4Z2
Protein Weight	23kDa
Shipping	Ice bag