

Anti-Neurog3 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y4Z2

Item Number: CM0022990

Introduction

NEUROG3 rabbit polyclonal antibody for WB and ELISA. Targets human neurogenin-3, a transcription factor critical for enteroendocrine and pancreatic endocrine cell differentiation. Ideal for developmental biology and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	NEUROG3 Rabbit pAb
Remarks/Alias	ngn3; Atoh5; NGN-3; Math4B; bHLHa7; G3
Species	Human
GeneID (Human)	50674
GeneID	50674
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human NEUROG3 (NP_066279.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y4Z2
Protein Weight	23kDa
Shipping	Ice bag