

## Anti-Neurogenin 2 (Ngn2) Monoclonal Antibody For Wb, Elisa - Q9H2A3

Item Number: CM0003547

### Introduction

Rabbit monoclonal antibody against Neurogenin 2, a neuronal transcription factor. Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for studying neuronal differentiation and development.

[Learn More](#)

| Parameter        | Value                                                                                                                        |
|------------------|------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Neurogenin 2 Rabbit mAb                                                                                                      |
| Remarks/Alias    | NGN2; Atoh4; ngn-2; Math4A; bHLHa8; Neurogenin 2                                                                             |
| Species          | Human                                                                                                                        |
| GeneID (Human)   | 63973                                                                                                                        |
| GeneID           | 63973                                                                                                                        |
| Immunogen        | Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Neurogenin 2 (Q9H2A3). |
| Source           | Rabbit                                                                                                                       |
| Category         | Monoclonal Antibodies                                                                                                        |
| Application      | WB, ELISA                                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                            |
| SWISS            | Q9H2A3                                                                                                                       |
| Protein Weight   | 29kDa                                                                                                                        |
| Shipping         | Ice bag                                                                                                                      |