

Anti-Neurogranin (Nrgn) Rabbit Polyclonal Antibody For Wb, Elisa - Q92686

Item Number: CM0027757

Introduction

Rabbit polyclonal antibody against human Neurogranin (NRGN), validated for WB and ELISA. Suitable for studies on synaptic remodeling, calmodulin binding, and Alzheimer disease. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	NRGN Rabbit pAb
Remarks/Alias	RC3; hng; NRGN
Species	Human
Gene ID (Human)	4900
Gene ID	4900
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-78 of human NRGN (NP_001119653.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92686
Protein Weight	8kDa
Shipping	Ice bag