

Anti-Neuronal Acetylcholine Receptor Subunit Alpha-4 (Acha4/Chrna4) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P43681

Item Number: CM0019753



Introduction

Rabbit polyclonal antibody against human ACHA4 (CHRNA4), validated for WB, IF/ICC, ELISA and cross-reacts with mouse. Ideal for studying neuronal nicotinic acetylcholine receptor function.

[Learn More](#)

Parameter	Value
Product Name	ACHA4 Rabbit pAb
Remarks/Alias	EBN; BFNC; EBN1; NACHR; NACRA4; NACHRA4; ACHA4
Species	Human
GenelD (Human)	1137
GenelD	1137
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-242 of human ACHA4 (P43681)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	P43681
Protein Weight	70kDa
Shipping	Ice bag