

Anti-Neuroplastin (Nptn) Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y639

Item Number: CM0026909



Introduction

Rabbit polyclonal antibody targeting human Neuroplastin (NPTN) validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studies on synaptic plasticity and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	NPTN Rabbit pAb
Remarks/Alias	GP55; GP65; SDR1; np55; np65; SDFR1; NPTN
Species	Human
Gene ID (Human)	27020
Gene ID	27020
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-200 of human NPTN (NP_059429.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y639
Protein Weight	44kDa
Shipping	Ice bag