


Anti-Neutrophil Elastase (Elane) Monoclonal Antibody For Wb, Ip, Elisa - P08246

Item Number: CM0009344



Introduction

High-quality rabbit anti-human Neutrophil Elastase (ELANE) monoclonal antibody optimized for WB, IP, and ELISA. A reliable tool for investigating neutrophil function and NETosis.

[Learn More](#)

Parameter	Value
Product Name	Neutrophil Elastase (ELANE) Rabbit mAb
Remarks/Alias	GE; NE; HLE; HNE; ELA2; SCN1; PMN-E; Neutrophil Elastase (ELANE)
Species	Human
GeneID (human)	1991
GeneID	1991
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 167-267 of human Neutrophil Elastase (ELANE) (P08246)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P08246
Protein Weight	29 kDa
Shipping	Ice bag