


Anti-Neutrophil Elastase (Elane) Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - P08246

Item Number: CM0016385



Introduction

Rabbit polyclonal antibody against human Neutrophil Elastase (ELANE) validated for WB, IF/ICC, IHC-P, ELISA; detects endogenous ELANE in human, mouse, rat samples; ideal for neutrophil biology, NETosis, and inflammatory disease research.

[Learn More](#)

Parameter	Value
Product Name	Neutrophil Elastase (ELANE) Rabbit pAb
Remarks/Alias	GE; NE; HLE; HNE; ELA2; SCN1; PMN-E; Neutrophil Elastase (ELANE)
Species	Human
GeneID (Human)	1991
GeneID	1991
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 30-267 of human Neutrophil Elastase (ELANE) (NP_001963.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08246
Protein Weight	29 kDa
Shipping	Ice bag