

Anti-Nexn Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q0Zgt2

Item Number: CM0016104

Introduction

Rabbit polyclonal antibody against human NEXN (Nexilin). Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for actin cytoskeleton and cardiac/skeletal muscle research.

[Learn More](#)

Parameter	Value
Product Name	NEXN Rabbit pAb
Remarks/Alias	CMH20; NELIN; NEXN
Species	Human
Gene ID (Human)	91624
Gene ID	91624
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 516-675 of human NEXN (NP_653174.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q0ZGT2
Protein Weight	81kDa
Shipping	Ice bag