


Anti-Nf-Kb P65/Rela Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa, Chip - Q04206

Item Number: CM0003400



Introduction

KO-validated rabbit monoclonal antibody targeting human NF-kB p65/RelA (Q04206). Suitable for WB, IHC-P, IF/ICC, ELISA, and ChIP applications with cross-reactivity to human, mouse, rat, and monkey.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NF-kB p65/RelA Rabbit mAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; NF-kB p65/RelA
Species	Human
GeneID (Human)	5970
GeneID	5970
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	Q04206
Protein Weight	60kDa
Shipping	Ice bag