

Anti-Nf-Kb P65/Rela Monoclonal Antibody (Ko Validated) For Wb, If/Icc, Ip, Elisa, Chip - Q04206

Item Number: CM0005449



Introduction

KO Validated rabbit monoclonal antibody targeting human NF-kB p65/RelA. Suitable for WB, IF/ICC, IP, ELISA, ChIP. Cross-reacts with human, mouse, rat, monkey. Ideal for inflammation, immunity, cancer research, and SARS-CoV-2 studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NF-kB p65/RelA Rabbit mAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; IA
Species	Human
GeneID (Human)	5970
GeneID	5970
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 450-551 of human NF-kB p65/RelA (NP_068810.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	Q04206
Protein Weight	58kDa/59kDa/60kDa
Shipping	Ice bag