

Anti-Nf-Kb P65/Rela Mouse Monoclonal Antibody For Wb, If/Icc, Elisa - Q04206

Item Number: CM0002114



Introduction

Mouse monoclonal antibody targeting human NF-kB p65/RelA (Q04206), validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for NF-kB signaling, inflammation, and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	NF-kB p65/RelA Mouse mAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; NF-kB p65/RelA
Species	Human
GeneID (Human)	5970
GeneID	5970
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 450-551 of human NF-kB p65/RelA (Q04206).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04206
Protein Weight	60kDa
Shipping	Ice bag