


Anti-Nf-Kb P65/Rela Rabbit Polyclonal Antibody For Wb - Q04206

Item Number: CM0016668



Introduction

Rabbit polyclonal antibody against human NF- κ B p65/RelA (SWISS: Q04206) validated for WB, IP, IF/ICC, and ELISA; cross-reacts with mouse and rat. Suitable for NF- κ B signaling and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	NF- κ B p65/RelA Rabbit pAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; NF- κ B p65/RelA
Species	Human
Gene ID (Human)	5970
Gene ID	5970
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 450-551 of human NF- κ B p65/RelA (NP_068810.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04206
Protein Weight	60 kDa
Shipping	Ice bag