


Anti-Nf-Kb2 Monoclonal Antibody For Wb, Elisa - Q00653

Item Number: CM0004541



Introduction

Rabbit monoclonal NF-κB2 antibody for WB and ELISA, recognizing human, mouse, and rat NF-κB2 (p52/p100). Involved in inflammation, immunity, and apoptosis. Ideal for NF-κB signaling research.

[Learn More](#)

Parameter	Value
Product Name	NF-κB2 Rabbit mAb
Remarks/Alias	p52; p100; H2TF1; LYT10; CVID10; LYT-10; NF-κB2; p49/p100; NF-κB2
Species	Human
GeneID (Human)	4791
GeneID	4791
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 800-900 of human NF-κB2 (NP_001070962.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00653
Protein Weight	97kDa
Shipping	Ice bag