


Anti-Nf-Kb2 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q00653

Item Number: CM0003763



Introduction

Rabbit monoclonal antibody against human NF-κB2 (p52/p100). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for studying NF-kappa-B signaling, inflammation, and immunity.

[Learn More](#)

Parameter	Value
Product Name	NF-κB2 Rabbit mAb
Remarks/Alias	p52; p100; H2TF1; LYT10; CVID10; LYT-10; NF-κB2; p49/p100; NF-κB2
Species	Human
GeneID (Human)	4791
GeneID	4791
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 745-877 of human NF-κB2 (Q00653)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00653
Protein Weight	97kDa
Shipping	Ice bag