


Anti-Nf-Kb2 (Phospho S866+S870) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q00653

Item Number: CM0001373



Introduction

Rabbit monoclonal phospho-specific antibody targeting NF-κB2 (p100/p52) at S866+S870, validated for WB, IHC-P, ELISA in human and mouse. Key tool for NF-κB signaling research, inflammation, and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NF-κB2-S866/S870 Rabbit mAb
Remarks/Alias	p52; p100; H2TF1; LYT10; CVID10; LYT-10; NF-κB2; p49/p100; Phospho-NF-κB2-S866/S870
Species	Human
GeneID (Human)	4791
Immunogen	Synthetic peptide around S866 & S870 of human NF-κB2
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q00653
Protein Weight	97kDa
Shipping	Ice bag