

Anti-Nfat2 Rabbit Monoclonal Antibody For Wb And Elisa - O95644

Item Number: CM0002109



Introduction

High-quality rabbit monoclonal antibody against human NFAT2 (NFATc1) for western blotting (WB) and ELISA applications. Suitable for research on T-cell signaling, cytokine expression, and osteoclastogenesis.

[Learn More](#)

Specification	Details
Product Name	NFAT2 Rabbit mAb
Remarks/Alias	NFAT2; NFATc; NF-ATC; NF-ATc1.2
Species	Human
GeneID (Human)	4772
GeneID	4772
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 844-943 of human NFAT2 (O95644).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O95644
Protein Weight	101kDa
Shipping	Ice bag