

Anti-Nfat2 Rabbit Monoclonal Antibody For Wb - 095644

Item Number: CM0004271



Introduction

Anti-NFAT2 Rabbit Monoclonal Antibody validated for WB and ELISA. Detects human and mouse NFAT2. Suitable for T-cell signaling, cytokine expression, and osteoclast differentiation research.

[Learn More](#)

Parameter	Specification
Product Name	NFAT2 Rabbit mAb
Remarks/Alias	NFAT2; NFATc; NF-ATC; NF-ATc1.2
Species	human
Gene ID (Human)	4772
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-305 of human NFAT2 (NP_001265598.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95644
Protein Weight	101kDa
Shipping	Ice bag