

Anti-Nfat2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O95644

Item Number: CM0019885



Introduction

NFAT2 rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. Recognizes human NFAT2 transcription factor; immunogen spans amino acids 115-298. Suitable for T-cell signaling and immunology research.

[Learn More](#)

Parameter	Value
Product Name	NFAT2 Rabbit pAb
Remarks/Alias	NFAT2; NFATc; NF-ATC; NF-ATc1.2
Species	Human
GeneID (Human)	4772
GeneID	4772
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 115-298 of NFAT2 (NP_001265598.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O95644
Protein Weight	101kDa
Shipping	Ice bag