


Anti-Nfatc2 Rabbit Polyclonal Antibody For Wb - Q13469

Item Number: CM0013980



Introduction

High-quality rabbit polyclonal antibody against human NFATC2, validated for WB and ELISA. Suitable for research on T-cell signaling, cytokine expression, and transcriptional regulation in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	NFATC2 Rabbit pAb
Remarks/Alias	JCOSL; NFAT1; NFATP; NFATC2
Species	Human
Gene ID (Human)	4773
Gene ID	4773
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13469
Protein Weight	100kDa
Shipping	Ice bag