

Anti-Nfatc3 Monoclonal Antibody For Wb, If-P, Elisa - Q12968

Item Number: CM0007269



Introduction

Rabbit monoclonal antibody against NFATC3 (NFAT4/NFATx), validated for Western blot, immunofluorescence, and ELISA. Detects human, mouse, and rat NFATC3 with a predicted molecular weight of 116 kDa. Ideal for T-cell signaling and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	NFATC3 Rabbit mAb
Remarks/Alias	NFAT4; NFATX; NF-AT4c; n339260; NFATC3
Species	Human
Gene ID	4775
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 850-1075 of human NFATC3 (NP_775188.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12968
Protein Weight	116kDa
Shipping	Ice bag