

Anti-Nfe2L1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14494

Item Number: CM0013909

Introduction

Rabbit polyclonal antibody against human NFE2L1 (NRF1/TCF11) for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat NFE2L1 (85 kDa). Ideal for oxidative stress, cholesterol metabolism, and proteostasis research.

[Learn More](#)

Parameter	Value
Product Name	NFE2L1 Rabbit pAb
Remarks/Alias	NRF1; NRF-1; TCF11; LCR-F1; NFE2L1
Species	Human
Gene ID (Human)	4779
Gene ID	4779
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 515-772 of human NFE2L1 (NP_003195.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14494
Protein Weight	85kDa
Shipping	Ice bag