

Anti-Nfe2L3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y4A8

Item Number: CM0015285

Introduction

Rabbit polyclonal antibody against human NFE2L3 (NRF3), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for erythroid-specific gene expression studies.

[Learn More](#)

Parameter	Value
Product Name	NFE2L3 Rabbit pAb
Remarks/Alias	NRF3; NFE2L3
Species	Human, Mouse, Rat
GeneID (human)	9603
GeneID	9603
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 330-550 of human NFE2L3 (NP_004280.5)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4A8
Protein Weight	76kDa
Shipping	Ice bag