

Anti-Nfia Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q12857

Item Number: CM0029907



Introduction

Rabbit polyclonal antibody against NFIA, a transcription factor and DNA replication regulator. Validated for WB, IHC-P, IF/ICC, and ELISA in human, mouse, and rat samples. For nuclear factor studies.

[Learn More](#)

Parameter	Value
Product Name	NFIA Rabbit pAb
Remarks/Alias	CTF; NF1-A; NFI-A; NFI-L; BRMUTD; NF-I/A; DEL1P32P31; C1DELp32p31; NFIA
Species	Human
GeneID (human)	4774
GeneID	4774
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 180-330 of human NFIA (NP_005586.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12857
Protein Weight	56kDa
Shipping	Ice bag