

Anti-Nfib/Nf1B2 Rabbit Monoclonal Antibody For Wb And Elisa - O00712

Item Number: CM0006145



Introduction

Rabbit monoclonal antibody against human NFIB/NF1B2, validated for WB and ELISA with cross-reactivity to human and mouse. Ideal for transcription regulation and brain development research.

[Learn More](#)

Parameter	Value
Product Name	NFIB / NF1B2 Rabbit mAb
Remarks/Alias	CTF; MACID; NF1-B; NFI-B; NFIB2; NFIB3; NF-I/B; NFI-RED; HMGIC/NFIB; NFIB / NF1B2
Species	Human
GeneID (human)	4781
GeneID	4781
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 321-420 of human NFIB / NF1B2 (NP_005587.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00712
Protein Weight	47kDa
Shipping	Ice bag