

Anti-Nfic Rabbit Polyclonal Antibody For Wb, Elisa - P08651

Item Number: CM0014757



Introduction

This rabbit polyclonal anti-NFIC antibody is validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Use for studying transcription regulation, DNA replication, and TERT gene activation.

[Learn More](#)

Parameter	Value
Product Name	NFIC Rabbit pAb
Remarks/Alias	CTF, NFI, CTF5, NF-I, NFIC
Species	Human
Gene ID	4782
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human NFIC (NP_005588.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08651
Protein Weight	56kDa
Shipping	Ice bag