


Anti-Nfil3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16649

Item Number: CM0013015



Introduction

High-quality anti-NFIL3 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Recognizes human and rat NFIL3 (Q16649), a circadian clock transcriptional regulator. Ideal for circadian rhythm and transcription research.

[Learn More](#)

Parameter	Value
Product Name	NFIL3 Rabbit pAb
Remarks/Alias	E4BP4; IL3BP1; NFIL3A; NF-IL3A; NFIL3
Species	Human
GeneID (Human)	4783
GeneID	4783
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 133-462 of human NFIL3 (NP_005375.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q16649
Protein Weight	51kDa
Shipping	Ice bag