

Anti-Nfkb1 Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa, Chip - P19838

Item Number: CM0008711



Introduction

Rabbit monoclonal NFKB1 antibody validated for WB, IHC-P, IP, ELISA, CHIP on human and mouse. Targets NF-kB p105/p50, a key transcription factor in immune and inflammatory responses. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	NFKB1 Rabbit mAb
Remarks/Alias	KBF1, EBP-1, NF-kB, CVID12, NF-kB1, NFKB-p50, Nfkapab, NF-kappaB, NFKB-p105, NF-kappa-B1, NF-kappabeta
Species	Human
Gene ID (Human)	4790
Gene ID	4790
Immunogen	Synthetic peptide derived from human NFKB1
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA, CHIP
Cross Reactivity	Human, Mouse
SWISS	P19838
Protein Weight	105 kDa
Shipping	Ice bag