

Anti-Nfkb1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P19838

Item Number: CM0014784



Introduction

Anti-NFKB1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human. Raised against recombinant NFKB1 (aa 740-964). Ideal for NF-kappa-B signaling and transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	NFKB1 Rabbit pAb
Remarks/Alias	KBF1; EBP-1; NF-kB; CVID12; NF-kB1; NFKB-p50; NFkappaB; NF-kappaB; NFKB-p105; NF-kappa-B1; NF-kappabeta; NFKB1
Species	Human
GeneID (human)	4790
GeneID	4790
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 740-964 of NFKB1 (NP_001158884.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	P19838
Protein Weight	105kDa
Shipping	Ice bag